


RoWalker

IMPROVE PATIENT MOBILITY
AND PREVENT FALLS

SPH
SAFE PATIENT HANDLING
MEDICAL



DISTRIBUTED BY:

eRGONOMIC
SAFETY
SOLUTIONS

For more information, please contact:

Matt Center Phone: (603) 748-9762
matt@ergonomicsafetyolutionsinc.com


www.ergonomicsafetyolutionsinc.com

March 2021

RoWalker

MOBILIZE PATIENTS EARLY AND SAFELY

- 400 lb weight capacity
- Mobilize Patients Early and Reduce Length of Stay
- Challenge patients to walk farther
- The RoWalker carries all the equipment needed for complex patients
- Reduce risk of falls



Twist knob below platform to widen or narrow platform width as desired



Release the hydraulic latch to allow platform adjustment. Close latch once platform is in desired position



Apply directional locks on both sides to keep device moving in straight path, apply directional lock to one side to correct a gait drift.

DISTRIBUTED BY:

For more information, please contact:

Matt Center Phone: (603) 748-9762
matt@ergonomicsafetyolutionsinc.com