

Epidural Positioning Device

SPH
SAFE PATIENT HANDLING
MEDICAL

The Epidural Positioning Device commonly known as the EPD, is a revolutionary device for positioning patients safely and securely during epidural procedures.

- **Reduces strain and injury risk to staff.**
- **Provides patient comfort with less invasion of their personal space compared to traditional positioning methods.**
- **Facilitates correct patient positioning by encouraging cervical, thoracic and lumbar flexion to maintain a solid and stable position.**

With the use of the positioning device, you will increase efficiency in operating rooms, labor and delivery suites, pain centers and any other facility/department where epidural and spinal anesthesia is performed.



Represented By:

e **ERGONOMIC**
SAFETY
SOLUTIONS

For more information, please contact:

Matt Center Phone: (603) 748-9762
matt@ergonomicsafety solutionsinc.com

www.ergonomicsafety solutionsinc.com

January 2019

180° Adjustable Face Rest

Tilts to attain and maintain proper cervical flexion



6 Position Adjustable Arm Rest

Adjustable range from the floor 32" - 43"
Allows patient to rest arms comfortably, provides added stability

Adjustable Torso Support

Increase lumbar and thoracic flexion with simple adjustments; chest support immobilizes back
(Pictured left: flat torso cushion
Pictured right: elevated torso cushion)

7 Position Adjustable Foot Rest

Adjustable range from the floor 10.5" - 19.5" accommodates a wide range of body sizes

High Quality Casters

Easily lock into place

**Sani-Cover™
Medical Grade
Disposable Face
Rest Covers**



For use with the E-EPD or any crescent-shaped face rest cushion.
The 50 gsm, spunbound fabric provides a comfortable surface for patients' faces and provides an aseptic barrier. The latex-free elastic allows for quick and easy changeover between patients. Available in a convenient, easy-dispenser box:
SPH-FRD-03-3110 - Box of 25 covers
SPH-FRD-03-3110-CASE - 8 Boxes of 25 ea.

Can be used on the EPD arm rests also.

EPD Specifications

Weight 12 lbs.
For Use With Patients Up To . . . 600 lbs.
Arm Rest Range 32" - 43"
Arm Rest Size 16.5"W x 10"D
Chest Piece Range 5.5" - 10"
Chest Piece Depth 3"
Chest Piece Size 8"W x 8"H
Shipping Weight 23 lbs.
Shipping Box . . . 33"L x 19"W x 15"H
Warranty Limited Lifetime

LDS Specifications (base)

Weight 42 lbs.
Foot Plate Size 20"W x 12"D
Foot Plate Range 10.5" - 19.5"
Width 20.5"
Length 23"
Shipping Weight 52 lbs.
Shipping Box . . . 30"L x 27"W x 21"H
Warranty Limited Lifetime

Represented By:



For more information, please contact:

Matt Center Phone: (603) 748-9762
matt@ergonomicsafety solutionsinc.com